


# PRELIMINARY ACTIVITY FOR Sugar Fermentation by Yeast

## Open Inquiry Version

Yeast can metabolize sugar in two ways, *aerobically*, with the aid of oxygen, or *anaerobically*, without oxygen. When yeast metabolizes a sugar under anaerobic conditions, ethanol ( $\text{CH}_3\text{CH}_2\text{OH}$ ) and carbon dioxide ( $\text{CO}_2$ ) gas are produced. An equation for the fermentation of the simple sugar glucose ( $\text{C}_6\text{H}_{12}\text{O}_6$ ) is:



The metabolic activity of yeast can be determined by the measurement of gas pressure inside the fermentation vessel.

In the Preliminary Activity, you will use a Gas Pressure Sensor to monitor the pressure inside a test tube as yeast metabolizes glucose anaerobically. When data collection is complete, you will perform a linear fit on the resultant graph to determine the fermentation rate.




Figure 1

After completing the Preliminary Activity, you will first use reference sources to find out more about sugar fermentation by yeast before you choose and investigate a researchable question dealing with fermentation. Some topics to consider in your reference search are:

- sugars
- glucose
- fermentation
- anaerobic respiration
- aerobic respiration
- yeast
- enzyme
- substrate
- enzyme inhibitor